



Fetal Monitor F-12L



Features

- 12.1inch high resolution touch screen, wide angel vieweing.Handwrien input of mother name, hospital name.
- Clinical event mark.
- Short time variance analysis..Probe of ffarlarming.
- Auto scoring and printout.
- Built-in rechargeable Li-Ion batery, support longgtime orking.
- 20000 hours data storage, support date checkup and review and printout.
- Clinical event analysis.
- WIFI networking with Luckcome central monitoring system
- Twins FHR overlap alarming.

Specifications

- FHR working frequency: 1.0MHZ
- Strength:<5 mW/cm²
- Transducer: Multi-cystals, wide wave beam, pulsed-wave working mode
- TOCO Range: 0-100 units
- Display:12.1" TFT touch screen
- Dimension:350D X 320D X85H (mm)
- Net Weight: 3.5kg
- Options:
 - Wins FHR measurement
 - Fetus waker
 - Support wired FHR/TOCO/FM monitoring
 - Support wired FHR/TOCO/FM/SPO2/NIBP/PR monitoringSupport wired FHR/TOCO/FM/SPO2/NIBP/PR/ECG/TEMP/RESP monitoring



ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)
 Ph.: +91-120-4166921, E-mail: info@acerbiomed.com
 Web.: www.acerbiomedicals.com, www.acerbiomed.com

Service Helpline
 18008894232

