

# Video Laryngoscope

## Model: VL-22



- Short Press to take photo/long press to record
- 3.5 inches touchable high definition screen.
- Ergonomic Design, Bacterio-static , Comfortable , and convient handle.
- Wireless Charging, compatible with various wireless chargers.
- High Light LED Supplementary Lighting Physical Heating anti- Fog technology.

The Very First Type-C Port, No distinction for up and side-down, reduce false operation

## Product Introduction

### Prospective Usage:

Clinically support to put patient's epiglottis up to expose glottis, precisely guide doctors to operate anesthesia through air tube or emergency treatment.

### Prospective Using Place:

Mainly can be used in department of anesthesiology, emergency department, ICU, and 120 urgent care center, etc.

### Product Advantage:

Small and portable device; high sensitive capacitance touchable screen, with multi-angle flexible; high clear digital technique, with integrated shooting and storage and function with physical heating and super antifog, helps the image to be clearly visible; the lithium battery has super long endurance and supports to wireless charging.

### Subtable Person:

Adult, Child, Neonatus



## Product Specification:

<b>Screen</b>	3.5 Inches LCD	<b>Field Angle</b>	> 60
<b>Rotational Angle</b>	Up & Down > 120 , Left & Right > 270	<b>Light Source</b>	LED Cold Light
<b>Capacity</b>	16G, it can shoot more than 100 thousand pictures or take video more than 20 hours (access to 256G)		
<b>Standby Time</b>	> 5 Hours		
<b>Charging Time</b>	< 3.5 Hours		