



Fetal Monitor F- 12 I



Features

- 12.1 inch high resolution touch screen, wide angle viewing. Handwritten input of mother name, hospital name.
- Clinical event mark.
- Short time variance analysis Probe of alarming.
- Auto scoring and printout.
- Built-in rechargeable Li-Ion battery, support long time working.
- 20000 hours data storage, support date checkup and review and printout.
- Clinical event analysis.
- Twins FHR overlap alarming.

Specifications

- FHR working frequency: 1.0MHZ
- Strength: <math><5 \text{ mW/cm}^2</math>
- Transducer: Multi- crystals, wide wave beam, pulsed-wave working mode
- TOCO Range: 0-100 units
- Display: 12.1" TFT touch screen
- Dimension: 350D X 320D X 85H (mm)
- Net Weight: 3.5kg
- Options: wins FHR measurement
Fetus waker
Support wireless FHR/TOCO/FM monitoring
Support wireless FHR/TOCO/FM/SPO2/NIBP/ PR/ECG/TEMP/RESP monitoring.



ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)
Ph.: +91-120-4166921, E-mail: info@acerbiomed.com
Web.: www.acerbiomedicals.com

Service Helpline
18008894232

